LED RL 4

LED-Ringlicht
LED Ring Light

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Art.-Nr. / Art.No.

G40 5750 000 / G40 5751 000 / G40 5752 000 / G40 5753 000 /
G40 5754 000 / G40 5755 000 / G40 5756 000 / G40 5757 000 /
G40 5758 000

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Vielen Dank

LED RL 4

This high quality product has been designed and developed to require a minimum of maintenance. Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

Warning

- Please disconnect plug from the main socket before cleaning the device.
- Storing and using the device only in dry rooms.

Safety Information

- To avoid eye damage, don't look directly into the light source when turned on.

Symbols

The following symbols are used in this manual.

CAUTION! This symbol indicates a potential risk and alerts you to proceed with caution.

Important: The evaluation of eye safety occurs according to the standard CIE S009/E:2002 („photobiological safety of lamps and lamp systems“). Within the risk grouping system of this CIE standard, the LEDs specified in this data sheet fall into the „low risk“ group (relating to devices in the visible spectrum with an exposure time of 100s). Under real circumstances (for exposure time, eye pupils, observation distance), it can be assumed that no endangerment to the eye exists from these devices. As a matter of principle, however, it should be mentioned that intense light sources have a high secondary exposure potential due to their blinding effect. As is also true when viewing other bright light sources (e.g. headlights), temporary reduction in visual acuity and afterimages can occur, leading to irritation, annoyance, visual impairment, and even accidents, depending on the situation.

Intended Use

The LED ring lights offer a standard diameter of 66mm. They have been mainly developed for lighting applications in the area of stereo microscopy. They provide a shadow free, constant and flicker free light and are basically suitable for all types of microscopes and cameras. With a ring adapter the diameter can be reduced in steps from 66mm down to 50 mm. The electronic light control allows a continuous and constant dimming from 0%-100%.

The ring light consists of 48 long life, pure-white or natural-white LEDs. Alternatively version in red, green, blue, IR and UV are offered. As an accessory a frosted glass diffuser creates a diffuse light. The LED ring lights with diffuser minimize mirror effects and can be used in quality control and inspection.

LED Ringlights RL-Series

LED Ring Light RL4	G40 5750 000
LED Ring Light RL4-NW	G40 5751 000
LED Ring Light RL4-R	G40 5752 000
LED Ring Light RL4-G	G40 5753 000
LED Ring Light RL4-B	G40 5754 000
LED Ring Light RL4-UV400	G40 5755 000
LED Ring Light RL4-UV375	G40 5756 000
LED Ring Light RL4-UV365	G40 5757 000
LED Ring Light RL4-IR	G40 5758 000

Package Contents

Description	Quantity
LED Ring Light RL4	1
Power Supply Plug, 24V with industry plug	1
Table Control	1

Getting Started

The LED ring light is a plug & play device very easy to use.

Please check before use whether ...

...the power cord is connected.

...the brightness is set to the desired value.

Description of Buttons and Switches

a) Electronic Dimming

For the continuous dimming of the light and for the complete switch off of the ring light.

b) Knurled Screws

To adapt the inner diameter of the ring light.

c) Operation Control

The blue LED glows as soon as the ring light has been switched on.


d) Connection from the ring light

Please connect correctly before use.

e) Connection from the power adaptor

Please connect correctly before use.

f) Power adaptor

LED RL 4

Troubleshooting

The information in the table below can help to resolve the possible operating problems quickly.

Problem	Causes	Solution
Power supply: No function at the device. No blue LED.	No network connection.	Please check, whether the power cord and all other cables are connected.

Storage, Maintenance, Repair

The device is maintenance-free. Some points, however you should take into consideration.

1. Please unplug the device before cleaning.
2. For the cleaning of the housing please use a smooth, semi-wet cloth, if necessary with a smooth household cleaning fluid.
3. No moisture should entering into the interior of the housing.
4. To avoid damages during the storing and transport please retain the packaging.
5. Repairs are only allowed by trained and experienced personnel. In case of improper handling the warranty expires.

Technical Service

Web Resources:

Visit our Qioptiq website at www.qioptiq.com for:

- Contact details of technical customer service.
- Qioptiq product overview and price list as well as infos about accessories.
- additional product information and special offers.

Contact:

For information or technical assistance contact directly Qioptiq Photonics GmbH & Co. KG.

Technical Specifications	
Weight:	382g
Inner diameter:	66mm
Outer diameter:	100mm
Operating Voltage:	100-240 VAC / 50-60 Hz
LED color:	white, neutral-white, green, blue, red, ultra violet, infra red
Color temperature:	typically 6.000-7.000 K, neutral-white: 4.600-5.600K
LED life span:	about 50.000h, UV-ring lights: about 10.000h (according to manufacturer information)
Light control:	continuous electronic dimming from 0%-100%
Material:	Aluminum anodised / plastic
Color of housing:	black
Working distance:	45mm-150mm

Accessories for RL4-Series		
	Frosted glass diffusor, inner-Ø 52mm	- Art.No.: G40 5764 000
	Polarizing Filters Set for RL 4	- Art.No.: G40 5762 000
	Reduction ring for adapting to Ø 50mm	- Art.No.: G40 5659 000
	Reduction ring for adapting to Ø 52mm	- Art.No.: G40 5660 000
	Reduction ring for adapting to Ø 54mm	- Art.No.: G40 5661 000
	Reduction ring for adapting to Ø 56mm	- Art.No.: G40 5662 000
	Reduction ring for adapting to Ø 58mm	- Art.No.: G40 5663 000
	Reduction ring for adapting to Ø 60mm	- Art.No.: G40 5664 000
	Reduction ring for adapting to Ø 62mm	- Art.No.: G40 5665 000
	Reduction ring for adapting to Ø 64mm	- Art.No.: G40 5666 000
	Protection glass	- Art.No.: G40 5765 000

LED RL 4

Warranty

Qioptiq Photonics GmbH & Co. KG warrants that this product will be free from defects in material and workmanship for a period of one (1) year from date of purchase. If a defect is present, Qioptiq will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products.

It is expressly agreed that this warranty will be in lieu of all warranties of fitness and in lieu of the warranty of merchantability.

Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead, it's your responsibility to correctly dispose of your equipment at lifecycle-end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

Notizen/Notes